


ALCONOX PRODUCT PRICE LIST

US SUGGESTED RETAIL PRICES EFFECTIVE JAN 1, 2018



Available Worldwide from leading laboratories, scientific, healthcare and industrial suppliers

Biodegradable. Concentrated. Interfering residue free. Alconox detergents set the standard for critical cleaning. Scientists, engineers and professionals specify Alconox detergents in their cleaning procedures for optimum economy, environmental safety and performance.



Alconox® Powdered Precision Cleaner

Original general-purpose powder for manual, soak, and ultrasonic cleaning of glass, metal, and plastic. No interfering residues. Has inhibitory residue test. USDA/FDA compliant. Dilute 1:100. pH 9.5

Size	Mfr Cat No	Suggested Retail US\$
4 lb box (1.8 kg)	1104-1	\$ 39.60
CS 9X4 lb (9X1.8 kg)	1104	\$ 297.00
25 lb box (11 kg)	1125	\$ 202.00
50 lb box (22.5 kg)	1150	\$ 377.00
100 lb drum (45 kg)	1101	\$ 678.00
300 lb drum (136 kg)	1103	\$ 1801.00
50X0.5 oz packs (50X15g)	1112-1	\$ 96.10
CS 12X50 packs	1112	\$ 887.00

Liquinox™ Critical Cleaning Liquid Detergent

Original general-purpose, phosphate-free liquid for manual, soak, and ultrasonic cleaning of glass, metal, and plastic. No interfering residues. Has inhibitory residue test. USDA/FDA compliant. Dilute 1:100. pH 8.5

1 qt bottle (0.95 L)	1232-1	\$ 23.30
CS 12X1 qt (12X.95 L)	1232	\$ 215.00
1 gal bottle (3.8 L)	1201-1	\$ 73.80
CS 4X1 gal (4X3.8 L)	1201	\$ 246.00
5 gal jerrycan (19 L)	1205	\$ 271.00
15 gal drum (57 L)	1215	\$ 724.00
55 gal drum (208 L)	1255	\$ 2320.00

Citranox® Liquid Acid Cleaner and Detergent

Liquid acid cleaner for manual, soak and ultrasonic cleaning or acid neutralizing. Corrosion-inhibited for glassware, metals, and plastic. Surfactant rinse aids for no interfering residues. Has inhibitory residue test. USDA and FDA compliant. Dilute 1:100. pH 2.5.

1 qt bottle (0.95 L)	1832-1	\$ 19.20
CS 12X1 qt (12X.95 L)	1832	\$ 192.00
1 gal bottle (3.8 L)	1801-1	\$ 60.90
CS 4X1 gal (4X3.8 L)	1801	\$ 203.00
5 gal jerrycan (19 L)	1805	\$ 254.00
15 gal drum (57 L)	1815	\$ 738.00
55 gal drum (208 L)	1855	\$ 2397.00

Detonox® Ultimate Precision Cleaner

For hand and ultrasonic use on difficult residues. Dilute 1:100. pH 10.5.

1 gal Bottle (3.8 L)	2301-1	\$ 88.10
CS 4x1 gal (4x3.8 L)	2301	\$ 271.00
5 gal jerrycan (19 L)	2305	\$ 301.00
55 gal drum (208 L)	2355	\$ 2702.00

Luminox® Low-Foaming Neutral pH Liquid Detergent

Low-foaming, neutral pH detergent for use in manual, ultrasonic, and machine cleaning. Safe-handling detergent used for noncorrosive, light-duty cleaning. No chelation or alkaline waste treat problems. Free rinsing, eliminating interfering residues. Dilute 3:100. pH 7.0

Size	Mfr Cat No	Suggested Retail US\$
1 gal bottle (3.8 L)	1901-1	\$ 84.90
CS 4X1 gal (4X3.8 L)	1901	\$ 283.00
5 gal jerrycan (19 L)	1905	\$ 358.00
15 gal drum (57 L)	1915	\$ 880.00
55 gal drum (208 L)	1955	\$ 2874.00

Tergazyme™ Enzyme-Active Powdered Detergent

Enzymatic cleaner for manual, soak, and ultrasonic cleaning of proteins, tissue, and body fluids from glass, metals, and plastic, with no interfering residues. Cleans filter membranes, fermenters and bioreactors. Has inhibitory residue test. FDA/USDA compliant. Dilute 1:100. pH 9.5

4 lb box (1.8 kg)	1304-1	\$ 42.30
CS 9X4 lb (9X1.8 kg)	1304	\$ 317.00
25 lb box (11 kg)	1325	\$ 200.00
50 lb box (22.5 kg)	1350	\$ 383.00
100 lb drum (45 kg)	1301	\$ 700.00
300 lb drum (136 kg)	1303	\$ 1843.00

Detergent 8® Low-Foaming Ion-Free Liquid Detergent

Nonionic detergent for PC board washing machines, manual, and ultrasonic cleaning. Ion- and chelate-free for no conductive residues on electronics and no chelate interference in waste. Phosphate free. Dilute 3:100. pH 11

1 gal Bottle (3.8 L)	1701-1	\$ 116.70
CS 4X1 gal (4X3.8 L)	1701	\$ 389.00
5 gal jerrycan (19 L) - Air OK	1705	\$ 502.00
15 gal drum (57 L)	1715	\$ 1205.00
55 gal drum (208 L)	1755	\$ 3910.00

Detojet™ Low-Foaming Liquid Detergent

Heavy-duty alkaline degreaser. Nonionic liquid for labware and machine washers, manual, and ultrasonic cleaning. Corrosion-inhibited for glass, metals, and plastic. No interfering residues. Has inhibitory residue test. USDA/FDA compliant. Dilute 1:200. pH 12

1 qt bottle (0.95 L)	1632-1	\$ 18.00
CS 12X1 qt (12X.95 L)	1632	\$ 166.00
1 gal bottle (3.8 L)	1601-1	\$ 56.10
CS 4X1 gal (4X3.8 L)	1601	\$ 187.00
5 gal can (19 L) - Air OK	1605	\$ 227.00
15 gal drum (57 L)	1615	\$ 627.00
30 gal drum (114 L)	1630	\$ 1131.00
55 gal drum (208 L)	1655	\$ 1793.00

Alcotabs® Tablet Pipet Detergent

Effervescent detergent tablet for maximum cleaning in siphon pipet washers. Slow release bubbles and free rinsing give no interfering residues. Has inhibitory residue test. Use 2 tabs/wash. pH 6.5

100 tablet bottle	1500-1	\$ 38.10
CS 6X100 bottles	1500	\$ 183.00

Alcojet® Low-Foaming Powdered Detergent

General-purpose nonionic powder for labware and machine washers, manual, and ultrasonic cleaning. Corrosion-inhibited for glassware, metals, and plastic. No interfering residues. Has inhibitory residue test. USDA authorized. Dilute 1:200. pH 11.5

Size	Mfr Cat No	Suggested Retail US\$
4 lb box (1.8 kg)	1404-1	\$ 35.60
CS 9X4 lb (9X1.8 kg)	1404	\$ 267.00
25 lb box (11 kg)	1425	\$ 179.00
50 lb box (22.5 kg)	1450	\$ 326.00
100 lb drum (45 kg)	1401	\$ 606.00
300 lb drum (136 kg)	1403	\$ 1663.00

Tergajet™ Low-Foaming Phosphate-Free Powder

Phosphate-free, nonionic powder for machine washers, hand, and ultrasonic cleaning. Corrosion-inhibited for labware. For water, microbiology, and environmental labs. No interfering residues. Has inhibitory residue test. USDA compliant. Dilute 1:200. pH 11.5

4 lb box (1.8 kg)	2204-1	\$ 39.50
CS 9X4 lb (9X1.8 kg)	2204	\$ 296.00
25 lb box (11 kg)	2225	\$ 192.00
50 lb box (22.5 kg)	2250	\$ 358.00
100 lb drum (45 kg)	2201	\$ 657.00
300 lb drum (136 kg)	2203	\$ 1765.00

Solujet™ Low-Foaming Phosphate-Free Liquid

Nonionic liquid for machine, spray-CIP, manual, soak, and ultrasonic cleaning. High emulsifying and penetrating formula. Corrosion-inhibited for glass, metals, and plastic. No interfering residues. Has inhibitory residue test. USDA and FDA compliant. Dilute 1:200. pH 12

1 gal bottle (3.8 L)	2101-1	\$ 66.60
CS 4X1 gal (4X3.8 L)	2101	\$ 205.00
5 gal can (19 L) - Air OK	2105	\$ 249.00
15 gal drum (57 L)	2115	\$ 703.00
30 gal drum (114 L)	2130	\$ 1278.00
55 gal drum (208 L)	2155	\$ 2214.00

Citrajel® Low-Foaming Liquid Acid Cleaner/Rinse

Liquid acid cleaner for machine, spray-CIP, manual, and ultrasonic cleaning or acid neutralizing. Corrosion-inhibited for glassware, metals, and plastic. Surfactant rinse for no interfering residues. Has inhibitory residue test. USDA/FDA compliant. Dilute 1:100. pH 2.5.

1 gal bottle (3.8 L)	2001-1	\$ 59.10
CS 4X1 gal (4X3.8 L)	2001	\$ 197.00
5 gal can (19 L)	2005	\$ 242.00
15 gal drum (57 L)	2015	\$ 676.00
30 gal drum (114 L)	2030	\$ 1291.00
55 gal drum (208 L)	2055	\$ 2243.00

Keylajet™ Low-Foaming High Alkaline Liquid

Corrosion inhibited formulation recommended for hard metals, steels, stainless steel, porcelain, ceramic, rigid plastics, cement, and elastomers. Free rinsing to give you reliable results and no interfering residues. Dilute 1:200. pH 13+

5 gal can (19 L) - Air OK	2405	\$ 265.00
55 gal drum (208 L)	2455	\$ 2447.00

Alconox, Inc.
Critical Cleaning Experts™

D E T E R G E N T S E L E C T I O N G U I D E

Application Key Concerns	Articles Cleaned/ Soil Removed	Cleaning Method	Recommended Cleaner Powder	Recommended Cleaner Liquid
Healthcare/Veterinary Effective preparation for sterilization, longer instrument life. Reduce waste.	Surgical, anaesthetic, and examining instruments and equipment. Catheters and tubes.	Manual, Ultrasonic, Soak	ALCONOX	LIQUINOX (p-free)
		Machine washer, sani-sterilizer	ALCOJET	DETOJET
	Blood, body fluids, tissue on instruments.	Manual, Ultrasonic, Soak	TERGAZYME	
Pharmaceutical/Medical Device/Biotechnology Passing cleaning validation for FDA good manufacturing practices. For stainless steel, glass, plastic, elastomer cleaning.	Titanium dioxide, petrolatum, oils, emulsions, ointments, carbopols, lacquers, zinc oxides, proteins, steroids, alcohols, sugars, and Eudragit* (L/S/L30/D55/NE30D) polymers.	Manual, Ultrasonic, Soak	ALCONOX	DETONOX
		Machine washer, power wash, CIP	ALCOJET TERGAJET (p-free)	KEYLAJET
	Inorganic residues, salts, metallics, pigments. Eudragit* (E/RL/RS/E100) polymers, amphoteric, coatings, amines, ethers, starches, alkaloids.	Manual, Ultrasonic, Soak		CITRANOX
		Machine washer, power wash, CIP		CITRAJET
	Protein/ferment residues. R/O, U/F membranes.	Manual, Ultrasonic, Soak	TERGAZYME	KEYLAJET
Laboratory/Environmental Reproducible results, no interfering residues, extending equipment life. Keep laboratory accreditation. Laboratory safety.	Glass, metal, plastic labware, ceramics, tissue culture, porcelain, clean rooms, animal cages, bioreactors, tubing, benches, safety equipment.	Manual, Ultrasonic, Soak	ALCONOX	LIQUINOX (p-free)
		Machine, power spray, labware washer, washer-sterilizer, cage-washers	ALCOJET TERGAJET (p-free)	DETOJET SOLUJET (p-free)
	Tubes, reusable pipets.	Siphon-type washer-rinsers	ALCOTABS (tablet)	
	Microbiology, water lab, and environmental sampling. Phosphate-sensitive labware. EPA procedures. (Acid for water rinse cycle.)	Field, Manual, Ultrasonic, Soak	TERGAJET (p-free)	LIQUINOX (p-free)
		Machine washer, labware washer	TERGAJET	SOLUJET-base CITRAJET-acid
	Radioactive equipment/contaminants.	Manual, Ultrasonic, Soak	ALCONOX	LIQUINOX
	Stopcock grease.	Machine washer, warewasher	ALCOJET	KEYLAJET
	Trace metals, metal oxides, scale, salts, starches, amines.	Manual, Ultrasonic, Soak		CITRANOX
		Machine washer, warewasher		CITRAJET
	Proteinaceous soils, bio-wastes, tissue, blood and other body fluids, fermentation residues.	Manual, Ultrasonic, Soak	TERGAZYME	
		Glassware washer	ALCOJET	SOLUJET
	Metalworking, Precision Manufacturing, and Optics Clean parts, avoid volatile solvents, strong acids and other hazardous chemicals.	Glass, ceramic, porcelain, stainless steel, plastic, rubber. Oils, chemicals, particulates.	Manual, Ultrasonic, Soak	ALCONOX
		Machine washer, power wash	ALCOJET TERGAJET (p-free)	DETOJET KEYLAJET
Aluminum, brass, copper, and other soft metal parts. Oils, chemicals, particulates (acid for oxides, salts, buffing compound).		Manual, Ultrasonic, Soak	ALCONOX TERGAJET (p-free)	LIQUINOX-base CITRANOX-acid
		Parts washer, power wash	ALCOJET TERGAJET (p-free)	SOLUJET-base CITRAJET-acid
Inorganics, metallic complexes, trace metals and oxides, scale, salts, metal brightening.		Manual, Ultrasonic, Soak		CITRANOX
		Parts washer, power wash		CITRAJET-acid KEYLAJET-base
Silicone oils, mold-release agents, buffing compounds.		Manual, Ultrasonic, Soak	ALCONOX	CITRANOX
		Parts washer, pressure spray	ALCOJET	KEYLAJET
	Delicate substrates/neutral for waste.	Manual, Ultrasonic, Soak Machine wash, pressure spray		LUMINOX (Neutral pH)
Electronics Avoid conductive residues, avoid CFCs, pass cleaning criteria.	Circuit boards, assemblies, screens, parts, conductive residues, resins, rosins, fluxes, particulates, salts.	Manual, Ultrasonic, Soak Machine washer, power spray board and screen washers		DETERGENT 8
	Ceramic insulators and components.	Manual, Ultrasonic, Soak	ALCONOX	LIQUINOX (p-free)
	Parts washers	ALCOJET	KEYLAJET	
Foodservice & Food Manufacturing Avoid interfering residues on food-contact equipment and effectively remove soils found in food operations.	Stainless steel, food production equipment.	Manual, Ultrasonic, Soak	ALCONOX	LIQUINOX
		Machine wash, pressure wash, CIP	ALCOJET	DETOJET
	Oxides, scale, trace metals, salts, milkstone and corrosion.	Manual, Ultrasonic, Soak		CITRANOX
		Machine wash, pressure wash, CIP		CITRAJET
	Filter membranes, proteins, biofouling.	Manual, Ultrasonic, Soak	TERGAZYME	KEYLAJET
Cosmetics Avoid cross-contamination.	Product contact surfaces (acids for pigments and salts).	Manual, Ultrasonic, Soak	ALCONOX	DETONOX CITRANOX-acid
		Parts washers, power spray	ALCOJET	KEYLAJET-base CITRAJET-acid
Nuclear Avoid waste interference	Reactor cavities, pipes, tools, protective equipment.	Manual, Soak, Spray		DETERGENT 8

p-free =phosphate-free

© 2018 Alconox, Inc.

* Eudragit is a ® registered trademark of Röhm GmbH & Co.

To purchase detergent contact your laboratory, healthcare, or industrial supply company.